
Potassium dihydrogen phosphate anhydrous, Hi-LR[®]

GRM1188

Product Identifier

CAS No.	:	7778-77-0
EC No.	:	231-913-4
Molecular Formula	:	KH ₂ PO ₄
Molecular Weight	:	136.09
Synonym	:	Potassium phosphate monobasic anhydrous, Hi-LR [®]
HS Code	:	2835 24 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to off-white crystals or powder or solid
Solubility	:	100 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	4.10 - 4.50
Chloride (Cl)	:	<= 0.01%
Sulphate (SO ₄)	:	<= 0.05%
Loss on drying (at 125 - 130°C, 2 hr)	:	<= 2.00%
Assay (NaOH Titration, on dry basis)	:	98.00 - 100.50%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	TC6615500
WGK	:	1