
EDTA ferric monosodium salt

GRM1197

Product Identifier

CAS No.	:	15708-41-5
EC No.	:	239-802-2
Molecular Formula	:	C ₁₀ H ₁₂ N ₂ NaFeO ₈
Molecular Weight	:	367.05
Synonym	:	Ethylenediaminetetraacetic acid ferric monosodium salt
HS Code	:	2922 49 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Light yellow to yellowish green to brown powder
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (Fe, Iodometry)	:	12.00 - 14.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AH4900000
WGK	:	2