

---

## Ammonium acetate, Hi-LR<sup>TM</sup>

**GRM1200**

### Product Identifier

---

CAS No.	:	631-61-8
EC No.	:	211-162-9
Molecular Formula	:	C <sub>2</sub> H <sub>7</sub> NO <sub>2</sub>
Molecular Weight	:	77.08
Synonym	:	Acetic acid ammonium salt, Hi-LR <sup>TM</sup>
HS Code	:	2915 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Colourless to off-white hygroscopic crystals or powder or solid
Solubility	:	1000 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	6.50 - 7.50
Chloride (Cl)	:	<= 0.01%
Iron (Fe)	:	<= 0.001%
Sulfate (SO <sub>4</sub> )	:	<= 0.02%
Assay (NaOH Back T/NT)	:	96.00 - 102.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AF3675000
WGK	:	1