

---

## Bees wax, white For histology

**GRM1213**

---

### Product Identifier

CAS No.	:	8012-89-3
EC No.	:	232-383-7
HS Code	:	1521 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to yellowish-white solid or pieces or plates, translucent whenthin, with a fine-grained, matt, non-crystalline fracture, when warmed becomes soft and pliable
Solubility	:	33.3 mg soluble in 1 mL of chloroform (Clear to hazy)
Acid value	:	5 - 15
Ester value	:	75 - 95
Saponification value	:	87 - 104

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1