
Bismuth (III) carbonate basic

GRM1221

Product Identifier

CAS No.	:	5892-10-4
EC No.	:	227-567-9
Molecular Formula	:	(BiO) ₂ CO ₃
Molecular Weight	:	509.97
HS Code	:	2836 99 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to faint beige powder
Solubility	:	Insoluble in water
Loss on drying (at 105°C, 2 hr)	:	<= 1.00%
Assay (Bi, EDTA Titration, on dry basis)	:	80.00 - 82.50%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3