

Boric acid, Hi-AR®

GRM1224


Product Identifier

CAS No.	:	10043-35-3
EC No.	:	233-139-2
EC Index No	:	005-007-00-2
Molecular Formula	:	H ₃ BO ₃
Molecular Weight	:	61.83
HS Code	:	2810 00 20
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white amorphous powder
Solubility	:	33.3 mg soluble in 1 mL of water
Chloride (Cl)	:	<= 0.001%
Calcium (Ca)	:	<=0.005%
Heavy metals (as Pb)	:	<= 0.001%
Iron (Fe)	:	<= 0.001%
Phosphate (PO ₄)	:	<= 0.001%
Sulfate (SO ₄)	:	<= 0.01%
Insoluble in methanol	:	<= 0.005%
Assay (NaOH Titration)	:	>= 99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H360FD
Precautionary Statement(s)	:	P201- P280- P308 + P313- P405- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	ED4550000
WGK	:	1

Other Information

Buffering agent