

GRM1226 Cadmium sulphate octahydrate**Product Number Packing**

GRM1226 : 500G

Product Information

Product Code : GRM1226
Product Name : Cadmium sulphate octahydrate
Molecular Formula : $3\text{CdSO}_4 \cdot 8\text{H}_2\text{O}$
Molecular Weight : 769.51
CAS No. : 7790-84-3
EC No. : 233-331-6
HS Code : 2833 29 20
EC Index No. : 048-009-00-9
Shelf Life : 4 years

Technical Specification

Appearance : White crystals or powder
Solubility : 33.3 mg soluble in 1 mL of water
pH (5% in water at 25°C) : 3.00 - 6.00
Chloride (Cl) : $\leq 0.005\%$
Iron (Fe) : $\leq 0.01\%$
Assay (EDTA Titration) : 98.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H301-H330-H340-H350-H360FD-H372-H410
Precautionary Statement(s) : P201-P260-P273-P284-P301+P310-P310
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 45-46-60-61-25-26-48/23/25-50/53
S-Phrase(s) : 53-45-60-61
WGK : 3
RTECS : EV2850000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 2570
Class : 6.1



Product Information

Packaging Group	:	3
Marine Pollutant	:	Yes
ADR/RID	:	2570 6.1/PG 3
IMDG	:	2570 6.1/PG 3
IATA	:	2570 6.1/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein