

---

## Glass wool

**GRM1232**

---

### Product Identifier

CAS No. : 65997-17-3  
EC No. : 266-046-0  
HS Code : 7019 80 00  
Storage : Below 30°C  
Shelf life : 4 years

---

### Technical Specification

Appearance : White coloured shining thin delicate fibres  
Solubility : Insoluble in water  
Alkalies (as Na<sub>2</sub>O) : ≤ 0.5%

---

### Safety Information

UN No. : Not dangerous goods  
Class : -  
Packing Group : -  
RTECS : LK3651000  
WGK : 1