

---

## Guar gum powder

**GRM1233**

---

### Product Identifier

CAS No.	:	9000-30-0
EC No.	:	232-536-8
HS Code	:	1302 32 90
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to beige to brown free flowing powder
Solubility	:	Dispersible in hot or cold water forming a colloidal solution
Starch	:	Passes test
Acid insoluble matter	:	<= 7.0%
Sulphated ash	:	<= 1.5%
Protein	:	<= 10.0%
Galactomannans	:	>= 70.00%
Loss on drying (at 105°C, 5 hr)	:	<= 15.0%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1