

GRM1234 Gum ghatti, Extra pure

Product Number

Packing

GRM1234	: 100G
GRM1234	: 500G

Product Information

Product Code	: GRM1234
Product Name	: Gum ghatti, Extra pure
CAS No.	: 9000-28-6
EC No.	: 232-534-7
HS Code	: 1301 90 00
Shelf Life	: 4 years

Technical Specification

Appearance	: Tan to dark brown lumps or granules
Solubility	: 1g is dispersed in 5 mL of water it forms a viscous, adhesive mucilage
Loss on drying (at 105°C, 5hr)	: <= 14.0%
Residue on ignition	: <= 6.0%

Risk and Safety Information

Storage Temperature(°C)	: Store below 30°C
-------------------------	--------------------

Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods