
SM Powder

GRM1254

Product Identifier

EC No.	:	310-127-6
HS Code	:	0402 10 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Pale yellow to cream amorphous homogeneous, free flowing powder
Solubility	:	33.3 mg soluble in 1 mL of water with slight opacity
pH of 2% w/v aqueous solution at 25°C	:	6.00 - 7.00
Water (K.F.)	:	<= 10.00%
Sulphated ash	:	<= 10.00%
Fat content	:	<= 1.50%
Total Nitrogen (anhydrous basis)	:	>= 4.70%
Total Protein (anhydrous basis)	:	>= 35.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1