

## Sodium benzoate, Hi-AR<sup>™</sup>

**GRM1260**


### Product Identifier

CAS No.	:	532-32-1
EC No.	:	208-534-8
Molecular Formula	:	C <sub>7</sub> H <sub>5</sub> NaO <sub>2</sub>
Molecular Weight	:	144.10
HS Code	:	2916 31 40
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White crystals or granules or powder or flakes
Solubility	:	33.3 mg soluble in 1 mL of water
Chloride (Cl)	:	<= 0.005%
Magnesium (Mg)	:	<= 0.005%
Calcium (Ca)	:	<= 0.01%
Heavy metals (as Pb)	:	<= 0.0001%
Loss on drying (at 105°C, 1 hr)	:	<= 1.0%
Iron (Fe)	:	<= 0.0001%
Sulfate (SO <sub>4</sub> )	:	<= 0.01%
Copper (Cu)	:	<= 0.0001%
Assay (NT, on dry basis)	:	min. 99.50

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319
Precautionary Statement(s)	:	P280- P305+P351+P338- P337 + P313
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DH6650000
WGK	:	1