
Sodium hexametaphosphate, Hi-AR®

GRM1282

Product Identifier

CAS No.	:	68915-31-1
EC No.	:	272-808-3
Synonym	:	Sodium polyphosphate, Hi-AR®
HS Code	:	2835 29 40
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white transparent crystals or platelets or granules or powder
Solubility	:	33.3 mg soluble in 1 mL of water
Chloride (Cl)	:	<= 0.03%
Iron (Fe)	:	<= 0.15%
Loss on ignition	:	<= 5.0%
Assay (as P ₂ O ₅ , NaOH Titration)	:	60.00 - 71.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	TR4950250
WGK	:	1