

---

## Propylparaben sodium salt

**GRM1290**

---

### Product Identifier

CAS No.	:	35285-69-9
EC No.	:	252-488-1
Molecular Formula	:	C <sub>10</sub> H <sub>11</sub> NaO <sub>3</sub>
Molecular Weight	:	202.18
Synonym	:	N-Propyl 4-hydroxy benzoate sodium salt
HS Code	:	2918 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to off-white hygroscopic crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
pH (0.1% in water at 25°C)	:	9.50 - 10.50
Chloride (Cl)	:	<= 0.033%
Sulfate (SO <sub>4</sub> )	:	<= 0.12%
Water (K.F.)	:	<= 5.00%
Assay (NT/Iodometry, on anhydrous basis)	:	99.00 - 102.00%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DH2850000
WGK	:	2