

## Methyl-p-hydroxybenzoate sodium salt

**GRM1291**


### Product Identifier

CAS No.	:	5026-62-0
EC No.	:	225-714-1
Molecular Formula	:	C <sub>8</sub> H <sub>7</sub> NaO <sub>3</sub>
Molecular Weight	:	174.13
Synonym	:	Methylparaben sodium salt
HS Code	:	2918 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White hygroscopic crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
pH (0.1% in water at 25°C)	:	9.50 - 10.50
Chloride (Cl)	:	<= 0.033%
Water (K.F.)	:	<= 5.0%
Sulfate (SO <sub>4</sub> )	:	<= 0.12%
Assay (NT/HPLC, on anhydrous basis)	:	95.00 - 102.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H315- H318- H412
Precautionary Statement(s)	:	P273- P280- P305+P351+P338- P310- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DH2455000
WGK	:	1