

2-Naphthol

GRM1292


Product Identifier

CAS No.	:	135-19-3
EC No.	:	205-182-7
Molecular Formula	:	C ₁₀ H ₈ O
Molecular Weight	:	144.17
Synonym	:	β-Naphthol
HS Code	:	2907 15 20
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or slightly pink crystals or powder or flakes or chips, turning brown on exposure to light
Solubility	:	100 mg soluble in 1 mL of ethanol
FTIR	:	Matches with the standard pattern
Melting range	:	120 - 125°C
Assay (GC)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302+H332- H400
Precautionary Statement(s)	:	P261- P270- P273- P301 + P312- P304 + P340- P501
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	QL2975000
WGK	:	2