



Product Information

GRM130

Egg yolk powder, From Hen

Product Number

Packing

GRM130 : 500G

Product Information

Product Code : GRM130
Product Name : Egg yolk powder, From Hen
Synonym : Pasteurised spray dried Hen Egg yolk powder
HS Code : 3507 90 90
Shelf Life : 4 years

Technical Specification

Appearance : Yellow powder
Solubility : Partially soluble in water
Protein : $\geq 30.0\%$
Water (K.F.) : $\leq 4.00\%$
Residue on ignition : $\leq 4.00\%$
Fat content : $\geq 57.00\%$
Mesh size : 100% passing through 18 mesh

Risk and Safety Information

WGK : 3
Storage Temperature(°C) : Store at 2 - 8°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein