

GRM1324 Mesalamine

Product Number Packing

GRM1324	: 100G
GRM1324	: 250G
GRM1324	: 50G

Product Information

Product Code	: GRM1324
Product Name	: Mesalamine
Synonym	: 5-Aminosalicylic acid
Molecular Formula	: C ₇ H ₇ N O ₃
Molecular Weight	: 153.14
CAS No.	: 89-57-6
EC No.	: 201-919-1
HS Code	: 2922 50 00
Shelf Life	: 4 years

Technical Specification

Appearance	: White to off-white to dark grey powder
Solubility	: 33.3 mg soluble in 1 mL of 1M hydrochloric acid
Melting range	: 275 - 280°C
FTIR (KBr disc)	: Matches with the standard pattern.
Assay (NaOH Titration)	: 95.00 - 101.00 %

GHS Safety Information

Hazard Statement(s)	: H315-H319-H335
Precautionary Statement(s)	: P261-P305+P351+P338
Signal Word	: Warning
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 36/37/38
S-Phrase(s)	: 26-36
WGK	: 2
RTECS	: VO1400000
Storage Temperature(°C)	: Store at 2 - 8°C
Hazard Symbol(s)	



Transport Information

Marine Pollutant	: -
------------------	-----



Product Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein