
Field's stain A powder, Hi-Cert™

GRM1328

Product Identifier

HS Code : 3804 19 29
Storage : Below 30°C
Shelf life : 4 years

Technical Specification

Appearance : Dark green to dark blue or black crystals or powder or solid
Solubility : 20 ppm in water yields clear solution
Absorption maxima 1 : 655 - 665 nm (20 ppm in water)
Absorption maxima 2 : 285 - 295 nm (20 ppm in water)

Safety Information

UN No. : Not dangerous goods
Class : -
Packing Group : -