

---

## Field's stain B powder, Hi-Cert™

**GRM1329**

---

### Product Identifier

HS Code : 3204 19 29  
Storage : Below 30°C  
Shelf life : 4 years

---

### Technical Specification

Appearance : Pinkish red powder or solid  
Solubility : 20 ppm in water yields clear solution  
Absorption maxima 1 : 515 - 525 nm (20 ppm in water)  
Absorption maxima 2 : 300 - 310 nm (20 ppm in water)  
Absorption maxima 3 : 250 - 260 nm (20 ppm in water)

---

### Safety Information

UN No. : Not dangerous goods  
Class : -  
Packing Group : -