

GRM1345 Benzamide**Product Number Packing**

GRM1345	: 100G
GRM1345	: 500G

Product Information

Product Code	: GRM1345
Product Name	: Benzamide
Synonym	: Benzoic acid amide
Molecular Formula	: C ₇ H ₇ NO
Molecular Weight	: 121.14
CAS No.	: 55-21-0
EC No.	: 200-227-7
HS Code	: 2924 29 98
Shelf Life	: 4 years

Technical Specification

Appearance	: White to off-white crystals or powder or chunks
Solubility	: 33.3 mg soluble in 1 mL of methanol
FTIR (KBr disc)	: Matches with the standard pattern
Melting range	: 124 - 130°C
Assay (GC/HPLC)	: min. 98.00%

GHS Safety Information

Hazard Statement(s)	: H302-H341
Precautionary Statement(s)	: P281
Signal Word	: Warning
Hazard Pictogram(s)	

**Risk and Safety Information**

R-Phrase(s)	: 22-68
S-Phrase(s)	: 36/37
WGK	: 1
RTECS	: CU8700000
Flash Point(°F)	: 356 °F
Flash Point(°C)	: 180°C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	





Product Information

Transport Information

Marine Pollutant	:	No
ADR/RID	:	Not Dangerous Goods
IMDG	:	Not Dangerous Goods
IATA	:	Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein