

Cadmium iodide

GRM1354


Product Identifier

CAS No.	:	7790-80-9
EC No.	:	232-223-6
EC Index No	:	048-007-00-8
Molecular Formula	:	CdI ₂
Molecular Weight	:	366.22
HS Code	:	2827 60 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to white to pale yellow or pearly white flakes or crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (AT/EDTA Titration)	:	90.00 - 101.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H301+H331- H351- H373- H410
Precautionary Statement(s)	:	P260- P273- P280- P301 + P310- P304 + P340- P403 + P233
UN No.	:	2570
Class	:	6.1
Packing Group	:	III
RTECS	:	EV1290000
WGK	:	3