

---

## Chitin, Flakes, From Shrimp shells

**GRM1356**

### Product Identifier

---

CAS No.	:	1398-61-4
EC No.	:	215-744-3
Molecular Formula	:	(C <sub>8</sub> H <sub>13</sub> NO <sub>5</sub> ) <sub>n</sub>
Synonym	:	Poly(N-acetyl-1,4-β-D-glucopyranosamine)
HS Code	:	3913 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White to brown flakes
Solubility	:	Insoluble in water
FTIR	:	Matches with the standard pattern.
Sulphated ash	:	<= 5.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1