

## Copper (I) iodide, Hi-LR™

**GRM1362**


### Product Identifier

CAS No.	:	7681-65-4
EC No.	:	231-674-6
Molecular Formula	:	CuI
Molecular Weight	:	190.45
Synonym	:	Cuprous iodide, Hi-LR™
HS Code	:	2827 60 10
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Off-white to grey-brown powder
Solubility	:	Insoluble in water
Assay (EDTA Titration)	:	99.00 - 102.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319- H335- H410
Precautionary Statement(s)	:	P261- P273- P305+P351+P338- P501
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	Not available
WGK	:	3