

Ammonium metavanadate

GRM1365


Product Identifier

| | | |
|-------------------|---|---------------------------------|
| CAS No. | : | 7803-55-6 |
| EC No. | : | 232-261-3 |
| Molecular Formula | : | NH ₄ VO ₃ |
| Molecular Weight | : | 116.98 |
| HS Code | : | 2841 90 00 |
| Storage | : | Below 30°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|----------------------------|---|---|
| Appearance | : | White to pale green or tan crystals or powder or solid |
| Solubility | : | 33.3 mg soluble in 1 mL of cooled mixture of water : sulphuric acid(150:50) |
| Chloride (Cl) | : | <= 0.2% |
| Sulfate (SO ₄) | : | <= 0.05% |
| Assay (FAS Titration) | : | min. 98.00% |

Safety Information

| | | |
|----------------------------|---|---|
| Hazard Pictogram(s) | : |  |
| Signal Word | : | Danger |
| Hazard Statement(s) | : | H301- H315- H319- H330- H335 |
| Precautionary Statement(s) | : | P260- P280- P284- P301 + P310- P304 + P340- P501 |
| UN No. | : | 2859 |
| Class | : | 6.1 |
| Packing Group | : | II |
| RTECS | : | YW0875000 |
| WGK | : | 3 |