

EDTA calcium disodium salt

GRM1371


Product Identifier

CAS No.	:	62-33-9
EC No.	:	200-529-9
Molecular Formula	:	C ₁₀ H ₁₂ N ₂ O ₈ CaNa ₂
Molecular Weight	:	374.27
Synonym	:	Ethylenediaminetetraacetic acid calcium disodium salt
HS Code	:	2922 49 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
pH (1% in water at 25°C)	:	6.50 - 8.00
Water (K.F.)	:	<= 13.00%
Assay (Bismuth nitrate T, on anhydrous basis)	:	92.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	EV7700000
WGK	:	3