

GRM1400 Potassium metabisulphite, Hi-AR™ (Potassium pyrosulphite, Hi-AR™)**Product Number Packing**

GRM1400 : 500G

Product Information

Product Code : GRM1400
Product Name : Potassium metabisulphite, Hi-AR™
Synonym : Potassium pyrosulphite, Hi-AR™
Molecular Formula : $K_2S_2O_5$
Molecular Weight : 222.32
CAS No. : 16731-55-8
EC No. : 240-795-3
Shelf Life : 4 years

Technical Specification

Appearance : White to light yellow crystals or powder
Solubility : 33.3 mg soluble in 1 mL of water
Arsenic (As) : $\leq 0.0001\%$
Chloride (Cl) : $\leq 0.005\%$
Copper (Cu) : $\leq 0.001\%$
Iron (Fe) : $\leq 0.001\%$
Lead (Pb) : $\leq 0.001\%$
Phosphate (PO₄) : $\leq 0.001\%$
Water insoluble matter : $\leq 0.005\%$
Zinc (Zn) : $\leq 0.001\%$
Assay (Iodometry) : 96.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H315-H319-H335
Precautionary Statement(s) : P261-P305 + P351 + P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 31-36/37/38
S-Phrase(s) : 26-36
WGK : 1
RTECS : TT4920000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)





Product Information

Transport Information

Marine Pollutant	:	No
ADR/RID	:	Not Dangerous Goods
IMDG	:	Not Dangerous Goods
IATA	:	Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein