
Potassium sulphate, Hi-AR™

GRM1403

Product Identifier

CAS No.	:	7778-80-5
EC No.	:	231-915-5
Molecular Formula	:	K ₂ SO ₄
Molecular Weight	:	174.26
HS Code	:	3104 30 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white to off-white crystals or powder or solid or granules
Solubility	:	33.3 mg soluble in 1 mL of water
Chloride (Cl)	:	<= 0.05%
Lead (Pb)	:	<= 0.003%
Iron (Fe)	:	<= 0.001%
Sodium (Na)	:	<= 1.0%
Assay (Cation exchange titration)	:	>= 98.50%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	TT5900000
WGK	:	1