

GRM1408 Silver chloride, Hi-LR™**Product Number Packing**

GRM1408 : 25G

Product Information

Product Code : GRM1408
Product Name : Silver chloride, Hi-LR™
Molecular Formula : AgCl
Molecular Weight : 143.32
CAS No. : 7783-90-6
EC No. : 232-033-3
Other Information : Extremely light sensitive
Shelf Life : 4 years

Technical Specification

Appearance : White to beige to grey powder or chunks.
Solubility : 33.3 mg soluble in 1 mL of ammonium solution
ICP : Confirms Silver component
Nitrate (NO₃) : ≤0.05%
Assay (Based on Trace Metal Analysis) : min. 99.00 %

GHS Safety Information

Hazard Statement(s) : H400
Precautionary Statement(s) : P273
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 50
S-Phrase(s) : 61
WGK : 3
RTECS : VW3563000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 3077
Class : 9
Packaging Group : 3



Product Information

Marine Pollutant	:	Yes
ADR/RID	:	3077 9/PG 3
IMDG	:	3077 9/PG 3
IATA	:	3077 9/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein