

GRM1413 Silver diethyldithiocarbamate, Hi-AR™/ACS**Product Number Packing**

GRM1413	: 25G
GRM1413	: 5G


Product Information

Product Code	: GRM1413
Product Name	: Silver diethyldithiocarbamate, Hi-AR™/ACS
Synonym	: Diethyldithiocarbamic acid silver salt, Hi-AR™/ACS
Molecular Formula	: C ₅ H ₁₀ AgNS ₂
Molecular Weight	: 256.14
CAS No.	: 1470-61-7
EC No.	: 216-003-7
HS Code	: 2843 29 00
Shelf Life	: 4 years


Technical Specification

Appearance	: White to yellow or green or khaki powder
Solubility	: 33.3 mg soluble in 1 mL of pyridine
FTIR (KBr disc)	: Matches with the standard pattern
Melting range	: 172 - 175°C
Assay (NH ₄ SCN Titration)	: 99.00 - 102.00%

GHS Safety Information

Hazard Statement(s)	: H315-H319-H335
Precautionary Statement(s)	: P261-P305+P351+P338
Signal Word	: Warning
Hazard Pictogram(s)	

Risk and Safety Information

R-Phrase(s)	: 36/37/38
S-Phrase(s)	: 26-36
WGK	: 3
Storage Temperature(°C)	: Store at 2 - 8°C
Hazard Symbol(s)	 Irritant

Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods



Product Information

IATA

: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein