
Silver metal powder

GRM1414

Product Identifier

CAS No.	:	7440-22-4
EC No.	:	231-131-3
Atomic Formula	:	Ag
Atomic Weight	:	107.87
HS Code	:	7106 10 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Silver to silver grey lustrous powder or granules
Solubility	:	33.3 mg soluble in 1 mL of hot nitric acid
Assay (NH ₄ SCN Titration)	:	min. 99.99%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	VW3500000
WGK	:	3