
Sodium thiosulphate anhydrous, Hi-AR™

GRM1420

Product Identifier

CAS No.	:	7772-98-7
EC No.	:	231-867-5
Molecular Formula	:	Na ₂ S ₂ O ₃
Molecular Weight	:	158.11
HS Code	:	2832 30 10
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
Loss on drying (at 105°C, 2hr)	:	<= 2.0%
Assay (Iodimetry, on dry basis)	:	97.00 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	XN6476000
WGK	:	1