
Sodium tungstate dihydrate, Hi-LR[™]

GRM1424

Product Identifier

CAS No.	:	10213-10-2
EC No.	:	236-743-4
Molecular Formula	:	Na ₂ WO ₄ ·2H ₂ O
Molecular Weight	:	329.85
HS Code	:	2841 80 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellow crystals or powder or flakes
Solubility	:	100 mg soluble in 1 mL of water
Chloride (Cl)	:	<= 0.01%
Iron (Fe)	:	<= 0.002%
Sulfate (SO ₄)	:	<= 0.02%
Assay (Lead nitrate Titration)	:	98.00 - 101.00%

Safety Information

Hazard Pictogram(s) :



Signal Word	:	Warning
Hazard Statement(s)	:	H302
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	YO7900000
WGK	:	1