

## DL-Tartaric acid, Hi-LR<sup>®</sup>

**GRM1429**


### Product Identifier

CAS No.	:	133-37-9
EC No.	:	205-105-7
Molecular Formula	:	C <sub>4</sub> H <sub>6</sub> O <sub>6</sub>
Molecular Weight	:	150.09
HS Code	:	2918 12 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White to off-white crystals or powder or solid
Solubility	:	100 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Specific rotation	:	-0.5° to +0.5° (c = 20% in water, at 20°C)
Assay (GC/NaOH Titration)	:	99.00 - 101.50%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H318
Precautionary Statement(s)	:	P280- P305+P351+P338- P310
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1