

GRM1444 p-Chloranil**Product Number Packing**

GRM1444 : 250G

Product Information

Product Code : GRM1444
Product Name : p-Chloranil
Synonym : Tetrachloro-1,4-benzoquinone; Tetrachloro-p-benzoquinone
Molecular Formula : $C_6Cl_4O_2$
Molecular Weight : 245.88
CAS No. : 118-75-2
EC No. : 204-274-4
HS Code : 2914 70 00
Shelf Life : 4 years

Technical Specification

Appearance : Yellow to yellow-green crystals or powder
Solubility : 33.3 mg soluble in 1 mL of hot toluene
FTIR (KBr disc) : Matches with the standard pattern
Melting range : ~289°C (dec.)
Assay (GC/Iodometry/HPLC) : 99.00 - 101.00%

GHS Safety Information

Hazard Statement(s) : H315-H319-H410
Precautionary Statement(s) : P273-P305+P351+P338-P501
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 36/38-50/53
S-Phrase(s) : 37-60-61
WGK : 2
RTECS : DK6825000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 3077
Class : 9
Packaging Group : 3



Product Information

Marine Pollutant	:	YES
ADR/RID	:	3077 9/PG 3
IMDG	:	3077 9/PG 3
IATA	:	3077 9/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein