

Phenazone

GRM1462


Product Identifier

CAS No.	:	60-80-0
EC No.	:	200-486-6
Molecular Formula	:	C ₁₁ H ₁₂ N ₂ O
Molecular Weight	:	188.23
Synonym	:	Antipyrine
HS Code	:	2933 11 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white to off-white crystals or powder
Solubility	:	1000 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Melting range	:	109 - 113°C
Assay (GC/NT)	:	99.00 - 101.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302
Precautionary Statement(s)	:	P264- P270- P301 + P312- P330- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	CD2450000
WGK	:	1