

---

## Tin (IV) oxide, Hi-AR™

**GRM1470**

---

### Product Identifier

CAS No.	:	18282-10-5
EC No.	:	242-159-0
Molecular Formula	:	SnO <sub>2</sub>
Molecular Weight	:	150.71
Synonym	:	Stannic oxide, Hi-AR™
HS Code	:	2825 90 10
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to light gray powder or chunks
Solubility	:	Insoluble in water
ICP	:	Confirms Tin component
Loss on ignition (at 1000°C, 1 hr)	:	<= 1.00%
Assay (Based on Trace Metal Analysis)	:	min. 99.00%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	XQ4000000
WGK	:	nwg