
Bilirubin, Hi-AR[™]

GRM1494

Product Identifier

CAS No.	:	635-65-4
EC No.	:	211-239-7
Molecular Formula	:	C ₃₃ H ₃₆ N ₄ O ₆
Molecular Weight	:	584.66
HS Code	:	2933 79 00
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Orange to brown crystals or powder
Solubility	:	10 ppm in chloroform yields clear to hazy yellow to orange solution
Absorption maxima	:	448 - 458 nm (10 ppm in chloroform)
Loss on drying (at 120°C, 2 hr)	:	<= 2.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DU3038000
WGK	:	2