

Hexamine, Hi-ARTM/ACS

GRM1499


Product Identifier

CAS No.	:	100-97-0
EC No.	:	202-905-8
EC Index No	:	612-101-00-2
Molecular Formula	:	C ₆ H ₁₂ N ₄
Molecular Weight	:	140.19
HS Code	:	2933 69 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder or solid
Solubility	:	100 mg soluble 1 mL of water
Heavy metals (as Pb)	:	<= 0.001%
Residue after ignition	:	<= 0.1%
Loss on drying (at 105°C, 4 hr)	:	<= 2.0%
Assay (H ₂ SO ₄ back-Titration, on dry basis)	:	min. 99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H228- H317
Precautionary Statement(s)	:	P210- P261- P280- P333 + P313- P362 + P364- P370 + P378
UN No.	:	1328
Class	:	4.1
Packing Group	:	III
RTECS	:	MN4725000
WGK	:	1