

---

## Silica gel, Self indicating blue coarse Water absorption capacity

**GRM151**

---

### Product Identifier

CAS No.	:	112926-00-8
EC No.	:	231-545-4
HS Code	:	2811 22 00
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	Blue solid or crystals or chunks
Solubility	:	Insoluble in water
Identification test	:	After absorbing moisture it becomes pink/colourless
Mesh Size	:	min. 90% retained on 6 mesh
Water adsorption capacity	:	>= 30.0% (24 hrs, 80% relative humidity)

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	nwg