
Molecular sieve 4Å, Hi-LR™

GRM1511


Product Identifier

CAS No.	:	70955-01-0
Molecular Formula	:	Na ₁₂ [(AlO ₂) ₁₂ (SiO ₂) ₁₂].xH ₂ O
Synonym	:	Sodium alumino silicate, Hi-LR™
HS Code	:	2842 10 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to brown beads
Solubility	:	Insoluble in water.

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3