

## Potassium chloroplatinate, ~40% Pt

**GRM1515**


### Product Identifier

CAS No.	:	16921-30-5
EC No.	:	240-979-3
Molecular Formula	:	$K_2PtCl_6$
Molecular Weight	:	485.99
Synonym	:	Potassium hexachloroplatinate (IV)
HS Code	:	2843 90 19
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Yellow to yellow-orange crystals or powder or chunks
Solubility	:	10 mg soluble in 10 mL of conc. nitric acid
Assay (Pt, Gravimetry)	:	39.40 - 40.60%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H301- H317- H318- H334
Precautionary Statement(s)	:	P261- P280- P301 + P310- P305+P351+P338- P342 + P311
UN No.	:	3290
Class	:	6.1
Packing Group	:	II
RTECS	:	TP1650000
WGK	:	1