

---

## 5% Palladium on activated carbon, Hi-LR<sup>®</sup>

**GRM1550**

---

### Product Identifier

CAS No.	:	7440-05-3
EC No.	:	231-115-6
Atomic Formula	:	Pd
Atomic Weight	:	106.42
HS Code	:	3815 12 10
Storage	:	Below 30°C
Shelf life	:	4 years


---

### Technical Specification

Appearance	:	Black powder
Solubility	:	Insoluble in water
Assay (Pd, Gravimetry)	:	4.80 - 5.20%

---

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H228
Precautionary Statement(s)	:	P210- P240- P241- P280- P370 + P378
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	nwg