
Palladium chloride

GRM1552


Product Identifier

CAS No.	:	7647-10-1
EC No.	:	231-596-2
Molecular Formula	:	PdCl ₂
Molecular Weight	:	177.33
HS Code	:	2843 90 19
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Orange-brown or red-brown crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of hydrochloric acid
Assay	:	min. 99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H301- H317
Precautionary Statement(s)	:	P280- P301 + P310
UN No.	:	3288
Class	:	6.1
Packing Group	:	III
RTECS	:	RT3500000
WGK	:	3