
Tween®40

GRM157

Product Identifier

CAS No.	:	9005-66-7
Synonym	:	Polyoxyethylenesorbitan monopalmitate
HS Code	:	3402 42 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Yellow to amber coloured viscous, oily liquid
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Hydroxyl value	:	89 - 105
Saponification value	:	41 - 52
Acid value	:	<= 2.0
Residue on ignition	:	<= 0.25%
Water (K.F.)	:	<= 3.0%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	WG2933000
WGK	:	1