

GRM1576

Irish moss

Product Number

Packing

GRM1576 : 100G
GRM1576 : 500G

Product Information

Product Code : GRM1576
Product Name : Irish moss
Synonym : I-Carrageenan
CAS No. : 9062-07-1
EC No. : 232-949-3
HS Code : 1302 39 00
Shelf Life : 4 years

Technical Specification

Appearance : White to light brown powder
Solubility : 10 mg soluble in 10 mL of hot water
pH (1% suspension at 25°C) : 8.00 - 11.00
Loss on drying (at 105°C, 2 hr) : <= 12.00%
Tota ash (on dry basis) : 15.00 - 40.00%
Viscosity : >= 5cP ((1.5% in water at 75°C, spindle No.1, 60 RPM)

Risk and Safety Information

WGK : 2
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods