
EDTA copper (II) disodium salt

GRM1680


Product Identifier

CAS No.	:	39208-15-6
EC No.	:	254-356-9
Molecular Formula	:	C ₁₀ H ₁₂ CuN ₂ Na ₂ O ₈
Molecular Weight	:	397.74
Synonym	:	Ethylenediaminetetraacetic acid copper (II) disodium salt
HS Code	:	2922 49 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Blue to dark blue to dark purple crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (Iodometry)	:	98.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1