
Poly(methyl methacrylate)

GRM1746

Product Identifier

| | | |
|-------------------|---|----------------------------------------------------------------|
| CAS No. | : | 9011-14-7 |
| Molecular Formula | : | $[\text{CH}_2\text{C}(\text{CH}_3)(\text{CO}_2\text{CH}_3)]_n$ |
| Molecular Weight | : | ~15,000 |
| HS Code | : | 3906 10 90 |
| Storage | : | Below 30°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|---------------------------------|---|--------------------------------------------------------------|
| Appearance | : | Colourless, transparent crystals or powder or lumps or beads |
| Solubility | : | 10 mg soluble in 1 mL of dichloromethane |
| Loss on drying (at 105°C, 2 hr) | : | <= 3.0% |

Safety Information

| | | |
|---------------|---|---------------------|
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | TR0400000 |
| WGK | : | NWG |