

Acetamide, Hi-AR™

GRM1780


Product Identifier

CAS No.	:	60-35-5
EC No.	:	200-473-5
EC Index No	:	616-022-00-4
Molecular Formula	:	C ₂ H ₅ NO
Molecular Weight	:	59.07
HS Code	:	2924 19 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white crystals or powder
Solubility	:	1000 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Melting range	:	78 - 81°C
Assay (GC)	:	min.99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H351
Precautionary Statement(s)	:	P201- P280- P308 + P313- P405- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AB4025000
WGK	:	1