
Soda lime with or without indicator, Hi-AR®/ACS

GRM1783

Product Identifier

CAS No.	:	8006-28-8
HS Code	:	3824 99 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to grey-white granules or pellets or chunks or solid
Description	:	Soda lime is a mixture of variable proportions of sodium hydroxide with calcium oxide or hydroxide. This reagent may or may not include an indicator
Solubility	:	Insoluble in water
Loss on drying (at 200°C, 18 hrs.)	:	<= 7.00%
Carbon dioxide absorption capacity	:	min. 19.00%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H314
Precautionary Statement(s)	:	P260- P280- P301+P330+P331- P303+P361+P353- P305+P351+P338- P501
UN No.	:	1907
Class	:	8
Packing Group	:	III
RTECS	:	Not available
WGK	:	1